




unyvero

Unyvero's sample-to-answer platform provides rapid results for severe infectious diseases in hospitalized patients

Powerful multiplex PCR technology combined with the broadest range of microorganism and resistance targets sets the Unyvero System apart.

The Unyvero System consists of:

- Lysator to lyse and process a variety of native samples
- Cockpit to manage testing process, display, store, and transmit results
- Analyzer to perform DNA testing with random-access, multiplex PCR



A single test handles one patient sample, analyzes over 100 DNA analytes and delivers reliable results within just 4-5 hours



Unyvero L4 Lysator



Unyvero C8 Cockpit



Unyvero A50 Analyzer



Unyvero is designed to expand with your growing needs

Applications for severe infections:

- Blood Culture – BCU
- Hospitalized Pneumonia – HPN
- Intra-Abdominal Infection – IAI
- Implant & Tissue Infection – ITI
- Urinary Tract Infection – UTI



unyvero

Curetis GmbH | Max-Eyth-Str. 42 | 71088 Holzgerlingen | Germany
Tel.: +49 (0)7031 49195 10 | Email: contact@curetis.com |

curetis



Blood Culture

Fast & Simple Syndromic Testing for Severe Infections - Improving Patient Outcomes

00587 Rev 4.0 © Curetis GmbH

Caution – Investigational Device. Limited by Federal US Law to Investigational Use. Not available for Sale in the United States.

www.curetis.com

> In industrialized countries, sepsis is responsible for as many deaths as heart attacks¹

- Sepsis is a major health issue and its recorded incidence is rising every year.^{2,3}
- Mortality rates can be as high as 50%.⁴
- Direct costs per sepsis patient range from €23,000 to €29,000 in Europe.⁵
- A microbiological analysis is recommended to narrow the empiric antimicrobial therapy as quickly as possible.⁶

Worldwide, someone dies of sepsis every 3-4 seconds.⁴

> Faster detection enables earlier optimization of therapy

The Unyvero BCU Application simultaneously identifies a large panel of bacteria, fungi and antibiotic resistance genes.

- Pathogen identification can take days using routine microbiology methods
- Every hour effective antibiotic treatment is delayed, sepsis mortality rate increases up to 8%⁷
- Early identification can help reduce morbidity and mortality, improve patient care and reduce healthcare

⁷ Kumar A *et al.*, Duration of hypotension before initiation of effective antimicrobial therapy is the critical determinant of survival in human septic shock. Crit Care Med. 2006; 34(6):1589-96

Unyvero Blood Culture (BCU) Cartridge

Gram-positive bacteria	Enterobacterales	Non-fermenting bacteria	Fungi	Resistance	Gene
<i>Staphylococcus aureus</i> <i>Coagulase negative staphylococci</i> <i>Streptococcus spp.</i> <i>Streptococcus agalactiae</i> <i>Streptococcus pneumoniae</i> <i>Streptococcus pyogenes/dysgalactiae</i> <i>Enterococcus spp.</i> <i>Enterococcus faecalis</i> <i>Listeria monocytogenes</i>	<i>Citrobacter freundii/koseri</i> <i>Escherichia coli</i> <i>Enterobacter cloacae complex</i> <i>Klebsiella aerogenes (E. aerogenes)</i> <i>Klebsiella oxytoca</i> <i>Klebsiella pneumoniae</i> <i>Klebsiella variicola</i> <i>Proteus spp.</i> <i>Serratia marcescens</i>	<i>Acinetobacter baumannii complex</i> <i>Pseudomonas aeruginosa</i> <i>Stenotrophomonas maltophilia</i>	<i>Aspergillus spp.</i> <i>Candida spp.</i> <i>Candida albicans</i> <i>Candida dubliniensis</i> <i>Candida glabrata</i> <i>I. orientalis (C. krusei)</i> <i>Candida parapsilosis</i> <i>Candida tropicalis</i>	Aminoglycoside Macrolide/Lincosamide Oxacillin Vancomycin 3rd generation Cephalosporins Carbapenem	<i>aac(6)/aph(2'')</i> <i>aacA4</i> <i>ermA</i> <i>mecA</i> <i>mecC</i> <i>vanA</i> <i>vanB</i> <i>ctx-M</i> <i>kpc</i> <i>imp</i> <i>ndm</i> <i>oxa-23</i> <i>oxa-24/40</i> <i>oxa-48</i> <i>oxa-58</i> <i>vim</i>
<i>Corynebacteriaceae</i>		Anaerobic bacteria			
<i>Corynebacterium spp.</i>		Cutibacterium acnes (P. acnes)			
		Other Gram-negative bacteria	<i>Mycobacteriaceae</i>		
		<i>Haemophilus influenzae</i> <i>Neisseria meningitidis</i>	<i>Mycobacterium spp.</i>		

> Clinical evidence demonstrates the benefits provided by the Unyvero solution

Study 1

Multicenter performance evaluation. Clinical laboratories from HDZ Bad Oeynhausen, UKE Hamburg and OWS Vienna.

Number of samples
178 positive blood cultures.

96.8% Sensitivity

99.8% Specificity

Value of resistance genes

- 119 resistance genes detected during study period.
- Important information for earlier targeted therapy and infection control.

Rapid identification is critical for survival

Using Unyvero, the average time to results was:

- Reduced by 11h compared to identification results.
- Reduced by 34h compared to full AST results.

Polymicrobial infections are severe

- 6/7 polymicrobial infections, with pathogens included in the panel, correctly identified.
- 5 samples with an additional micro-organism were detected using Unyvero.

11 hours saved



34 hours saved



Burrack-Lange *et al.*, Multicenter assessment of the rapid Unyvero Blood Culture molecular assay J. Med. Microbiology 2018; 67(9): 1294-1301
ID: Identification
AST: Antibiotic Susceptibility Testing

Study 2

Comparison with routine microbiology. University Hospital Essen, Germany.

Number of samples
50 positive blood cultures.
7 spiked blood cultures.

Study population

- 46 patients.
- Age 1 month-45 years.

100% Agreement for resistance markers

Spiked cultures

- 7/7 Spiked bacteria correctly identified.
- 3/3 Resistance markers correctly identified (*mecC*, *vanA* and *vanB*).



mecC



vanA



vanB

Conclusion

The Unyvero BCU Application is a useful tool for the rapid detection of pathogens and resistance markers directly from positive blood cultures.

Schmidt *et al.*, (2018) poster presentation DGHM.



Easy Workflow



Multiple Sample Types



24/7 Results